

Clavicle Fixation System

INSTRUCTIONS FOR USE FOR THE UNITED STATES OF AMERICA

R: For use by physicians only. Federal Law restricts this device to sale by or on the order of a physician.

This package insert is designed to provide Instructions for Use of the Clavicle Fixation System. **Failure to follow instructions may lead to patient injury.**

Description

The Clavicle Fixation System features bone plates for the repair of the superior midshaft and superior lateral regions, as well as the anterior midshaft and anterior lateral areas. This system can also facilitate a dual plate repair. Included in the Clavicle Fixation system are titanium (Ti) bone screws and pegs, Cobalt-Chrome (CoCr) bone screws, K-wires, and specialized instrumentation. The Clavicle Fixation System plates are available in various sizes and are made of medical-grade titanium alloy. The system is provided non-sterile and is sterilized at the user's facility.

The Clavicle Fixation System is comprised of:

- Titanium alloy plates, pegs, and screws (CoCr screws also available)
- Stainless steel K-wires (for provisional fixation; not for implantation)
- System-specific instrumentation

Skeletal Dynamics Clavicle Fixation Plates do not incorporate a medicinal substance.

The implants and associated instrumentation are provided non-sterile and are intended for steam sterilization at the user facility.

Indications

The Skeletal Dynamics Clavicle Fixation System is indicated for fixation of fractures, osteotomies and non-unions of the clavicle including osteopenic bone.

Contraindications

Prior to using the Clavicle Fixation System, ensure that none of the following patient conditions are present: active or latent infection, insufficient quantity or quality of bone and/or soft tissue, material sensitivity, or patients who are unwilling or incapable of following post operative care instructions. The use of the Clavicle Fixation System is contraindicated in cases where the implant(s) would cross open epiphyseal plates in skeletally immature patients.

⚠ Warnings & Precautions

- For safe effective use of the implant, the surgeon must be thoroughly familiar with the surgical technique, implant, and associated instruments. Careful consideration is required to avoid injury to the acromioclavicular, coracoclavicular, and sternoclavicular ligaments, as soft-tissue damage may already be present.
- The lungs, subclavian artery and brachial plexus are located in close proximity to the clavicle and are at risk during drilling. Use extreme caution, and when necessary, place a Crego Elevator beneath the clavicle to protect against drill penetration beyond the far cortex.
- Use caution when drilling with K-Wires or drills as it may result in vascular or nerve injury.
- When drilling in an oblong hole to adjust plate position, improper sequencing or technique may compromise fracture reduction. Ensure correct placement of screws in round or oblong holes based on the desired compression across the fracture.

- To ensure proper screw trajectory, the fixed-angle drill guide must be used when screw placement is required along the axis of a threaded hole. Failure to do so may result in misalignment or inadequate fixation.
- If the screws are not attached and/or fully tightened, a non-union, delayed union or construct failure may occur.
- The Post components MUST be used with only one 2.5mm Variable Angle Locking (VAL) Big Head screw. Use of any other screw diameter or screw type with the post is contraindicated. Use of incompatible screws may result in improper fit, unintended contact between components, fracture site distraction, increased mechanical stress, loss of fixation, or device failure.
- The patient should be informed about the importance of following the postoperative rehabilitation prescribed to fully understand the limitations in activities of daily living. The patient must be warned that failure to follow postoperative care instructions may cause the implant or treatment to fail.
- The patient must be cautioned about the use, limitations, and potential adverse effects of this device including the possibility of delayed union, non-union, device or treatment failure as a result of loosening, stress, excessive activity, or weight bearing or load bearing, and the possibility of nerve or soft tissue damage related to either surgical trauma or the presence of the device.
- Potential Clavicle Fixation System construct failures, such as implant breakage, loosening, fixation failure, delayed union, or non-union, may occur as a result of noncompliance with postoperative rehabilitation, excessive activities, or construct overloading.
- DO NOT reuse any of the Clavicle Fixation System implantable components. Reuse may compromise the structural integrity of the construct.
- Protect the Clavicle Fixation System's implantable components against scratching or nicking. Such stress concentration can lead to implant failure.
- Before using the Clavicle Fixation System, inspect all implants and instruments for wear, disfiguration, and physical damage. If evidence of wear, disfiguration or physical damage is found, DO NOT use and contact your local Skeletal Dynamics representative or the Skeletal Dynamics Customer Care Department.
- DO NOT permanently implant the Skeletal Dynamics K-Wires; they are intended to be used during provisional fixation of the Clavicle Plating. Additionally, improper placement may interfere with final fixation or compromise implant stability.
- DO NOT use screw lengths that will excessively protrude through the far cortex, as it may result in soft tissue irritation.
- Caution should be taken when contouring plates. Bending the plates may weaken or break the plates.
- The maximum angulation of the variable angle locking screw should not exceed 10° from the trajectory of the respective hole.
- The Clavicle Fixation System is to be used only with Skeletal Dynamics instruments, implants, and accessories.
- Dispose of contaminated implants and instruments per established facility guidelines and protocols.
- Caution should be taken for interference to pacemakers during electrocautery or by uncertified drill power sources.
- Seek medical help immediately if the implant malfunctions.
- Excess angulation may compromise plate integrity or screw purchase.
- To maintain traceability of the Clavicle Fixation System implantable components, you must record each of the respective components' LOT numbers into the patient's records.

Potential Adverse Events

The following are potential risks that have been associated with orthopedic surgery: infection, nonunion, persistent pain, damage to nerves, bone, blood vessels, or soft tissue, abnormal bone growth, discomfort or abnormal sensations, reduced range of motion, implant breakage, loosening or migration of the implants resulting in misalignment.

⚠ MR Safety Information

The devices included in the Clavicle Plating System have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Clavicle Plating System in the MR environment is unknown. Scanning a patient who has this device implanted may result in patient injury.

Directions for Use

The Clavicle Fixation System should only be used by surgeons who have experience with this system. Each surgeon must evaluate the appropriateness for the use of the Clavicle Fixation System based on their clinical experiences. Please refer to the Clavicle Fixation System's Surgical Technique Guide (MKT-00478-00) to review the surgical approach as described by Jorge L. Orbay, M.D. of the *Miami Hand and Upper Extremity Institute* located in Miami, Florida.

Cleaning, Sterilization, and Inspection

For instructions on cleaning, disinfection, sterilization, and inspection of the Clavicle Fixation System products please refer to the cleaning and sterilization instructions for use (IFU-04056-00). Please refer to the Clavicle Fixation System Surgical Technique Guide (MKT-00478-00) for proper tray arrangement.

Resources

The latest version of the Instructions for Use may be requested in a physical format by email (orders@skeletaldynamics.com) or by phone (+1-877-753-5396). The physical copy will be provided within 7 calendar days of receiving a request from the user or at the time of delivery of the device, if so, requested at the time of order.

For the most current instructions for use visit www.skeletaldynamics.com/resources. Instructions for Use should always be reviewed before using or implanting a device.

Disclaimer of Warranty and Limited Remedies

Skeletal Dynamics, Inc makes no express or implied warranty, including any implied warranty of merchantability or fitness for a particular purpose, on the product(s) described in this publication. Skeletal Dynamics, Inc. shall not be liable under any circumstances for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has authority to bind Skeletal Dynamics, Inc. to any representation or warranty except as specifically System forth in this publication. Descriptions or specifications provided by Skeletal Dynamics, Inc. in any publication is only included to generally describe the product when manufactured and do not constitute any express warranties.



Skeletal Dynamics, Inc.
7300 N. Kendall Dr. / Suite 800
United States of America,
Miami, FL 33156
1-877-753-5396



SYMBOLS GLOSSARY

Symbols that follow BS EN ISO 15223-1 / Medical Devices-Symbols to be used with medical device labels, labelling and information to be supplied.

Symbol	Symbol Reference Number and Title	Description
	5.2.7 Non-sterile	Indicates a medical device that has not been subjected to a sterilization process
	5.2.6 Do not resterilize	Indicates a medical device that is not to be resterilized
	5.1.6 Catalogue Number	Indicates the manufacturer's catalogue number so the medical device can be identified
	5.4.3 Consult Instructions for Use	Indicates the need for the user to consult the instructions for use
	5.1.5 Batch Code	Indicates the manufacturer's batch code or lot can be identified
	5.4.2 Do Not Re-use	Indicates a medical device that is intended for one single use only
	5.1.1 Manufacturer	Indicates the medical device manufacturer
	5.2.8 Do Not Use If Package is Damaged	Indicates a medical device that should not be used if the package has been damaged or opened
	5.1.4 Use-by-date	Indicates the date after which the medical device is not to be used
	5.7.10 Unique Device Identifier	Indicates a carrier that contains unique device identifier information
	5.7.7 Medical device	Indicated the item is a medical device
	5.1.3 Date of Manufacture	Indicates the date when the medical device was manufactured

Clavicle Fixation System Tray

Inventory Control Sheet

Titanium Alloy Superior Midshaft Clavicle Plates

Clavicle Plate, Superior Midshaft, 3.5mm, Size 1, 8 Hole, Right, Ti CPSM-35-8S1R (01)00841506142838  (01) 00841506142838	Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 13 Hole, Right, Ti CPSM-35-13S5R (01)00841506142975  (01) 00841506142975
Clavicle Plate, Superior Midshaft, 3.5mm, Size 1, 8 Hole, Left, Ti CPSM-35-8S1L (01)00841506142845  (01) 00841506142845	Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 13 Hole, Left, Ti CPSM-35-13S5L (01)00841506142982  (01) 00841506142982
Clavicle Plate, Superior Midshaft, 3.5mm, Size 1, 9 Hole, Right, Ti CPSM-35-9S1R (01)00841506147550  (01) 00841506147550	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 8 Hole, Right, Ti CPSM-27-8S1R (01)00841506142999  (01) 00841506142999
Clavicle Plate, Superior Midshaft, 3.5mm, Size 1, 9 Hole, Left, Ti CPSM-35-9S1L (01)00841506147567  (01) 00841506147567	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 8 Hole, Left, Ti CPSM-27-8S1L (01)00841506143002  (01) 00841506143002
Clavicle Plate, Superior Midshaft, 3.5mm, Size 2, 9 Hole, Right, Ti CPSM-35-9S2R (01)00841506142852  (01) 00841506142852	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 9 Hole, Right, Ti CPSM-27-9S1R (01)00841506143019  (01) 00841506143019
Clavicle Plate, Superior Midshaft, 3.5mm, Size 2, 9 Hole, Left, Ti CPSM-35-9S2L (01)00841506142869  (01) 00841506142869	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 9 Hole, Left, Ti CPSM-27-9S1L (01)00841506143026  (01) 00841506143026
Clavicle Plate, Superior Midshaft, 3.5mm, Size 2, 10 Hole, Right, Ti CPSM-35-10S2R (01)00841506147574  (01) 00841506147574	Clavicle Plate, Superior Midshaft, 2.7mm, Size 2, 9 Hole, Right, Ti CPSM-27-9S2R (01)00841506143033  (01) 00841506143033
Clavicle Plate, Superior Midshaft, 3.5mm, Size 2, 10 Hole, Left, Ti CPSM-35-10S2L (01)00841506147581  (01) 00841506147581	Clavicle Plate, Superior Midshaft, 2.7mm, Size 2, 9 Hole, Left, Ti CPSM-27-9S2L (01)00841506143040  (01) 00841506143040
Clavicle Plate, Superior Midshaft, 3.5mm, Size 3, 10 Hole, Right, Ti CPSM-35-10S3R (01)00841506142876  (01) 00841506142876	Clavicle Plate, Superior Midshaft, 2.7mm, Size 2, 10 Hole, Right, Ti CPSM-27-10S2R (01)00841506143057  (01) 00841506143057
Clavicle Plate, Superior Midshaft, 3.5mm, Size 3, 10 Hole, Left, Ti CPSM-35-10S3L (01)00841506142883  (01) 00841506142883	Clavicle Plate, Superior Midshaft, 2.7mm, Size 2, 10 Hole, Left, Ti CPSM-27-10S2L (01)00841506143064  (01) 00841506143064
Clavicle Plate, Superior Midshaft, 3.5mm, Size 3, 11 Hole, Right, Ti CPSM-35-11S3R (01)00841506142890  (01) 00841506142890	Clavicle Plate, Superior Midshaft, 2.7mm, Size 3, 10 Hole, Right, Ti CPSM-27-10S3R (01)00841506143071  (01) 00841506143071
Clavicle Plate, Superior Midshaft, 3.5mm, Size 3, 11 Hole, Left, Ti CPSM-35-11S3L (01)00841506142906  (01) 00841506142906	Clavicle Plate, Superior Midshaft, 2.7mm, Size 3, 10 Hole, Left, Ti CPSM-27-10S3L (01)00841506143088  (01) 00841506143088
Clavicle Plate, Superior Midshaft, 3.5mm, Size 4, 11 Hole, Right, Ti CPSM-35-11S4R (01)00841506142913  (01) 00841506142913	Clavicle Plate, Superior Midshaft, 2.7mm, Size 3, 11 Hole, Right, Ti CPSM-27-11S3R (01)00841506143095  (01) 00841506143095

Clavicle Fixation System Tray

Inventory Control Sheet

Clavicle Plate, Superior Midshaft, 3.5mm, Size 4, 11 Hole, Left, Ti CPSM-35-11S4L (01)00841506142920  (01) 00841506142920	Clavicle Plate, Superior Midshaft, 2.7mm, Size 3, 11 Hole, Left, Ti CPSM-27-11S3L (01)00841506143101  (01) 00841506143101
Clavicle Plate, Superior Midshaft, 3.5mm, Size 4, 12 Hole, Right, Ti CPSM-35-12S4R (01)00841506142937  (01) 00841506142937	Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 13 Hole, Right, Ti CPSM-35-13S5R (01)00841506142975  (01) 00841506142975
Clavicle Plate, Superior Midshaft, 3.5mm, Size 4, 12 Hole, Left, Ti CPSM-35-12S4L (01)00841506142944  (01) 00841506142944	Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 13 Hole, Left, Ti CPSM-35-13S5L (01)00841506142982  (01) 00841506142982
Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 12 Hole, Right, Ti CPSM-35-12S5R (01)00841506142951  (01) 00841506142951	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 8 Hole, Right, Ti CPSM-27-8S1R (01)00841506142999  (01) 00841506142999
Clavicle Plate, Superior Midshaft, 3.5mm, Size 5, 12 Hole, Left, Ti CPSM-35-12S5L (01)00841506142968  (01) 00841506142968	Clavicle Plate, Superior Midshaft, 2.7mm, Size 1, 8 Hole, Left, Ti CPSM-27-8S1L (01)00841506143002  (01) 00841506143002
Titanium Alloy Superior Cluster Clavicle Plates	
Clavicle Plate, Superior Cluster, 2.7mm, 12 Hole, Right, Ti CPSC-27-12R (01)00841506143118  (01) 00841506143118	Clavicle Plate, Superior Cluster, 3.5mm, 12 Hole, Right, Ti CPSC-35-12R (01)00841506143132  (01) 00841506143132
Clavicle Plate, Superior Cluster, 2.7mm, 12 Hole, Left, Ti CPSC-27-12L (01)00841506143125  (01) 00841506143125	Clavicle Plate, Superior Cluster, 3.5mm, 12 Hole, Left, Ti CPSC-35-12L (01)00841506143149  (01) 00841506143149
Titanium Alloy Superior Lateral Clavicle Plates	
Clavicle Plate, Superior Lateral, 2.7mm, 7 Hole, Right, Ti CPSL-27-7R (01)00841506143156  (01) 00841506143156	Clavicle Plate, Superior Lateral, 3.5mm, 8 Hole, Right, Ti CPSL-35-8R (01)00841506143170  (01) 00841506143170
Clavicle Plate, Superior Lateral, 2.7mm, 7 Hole, Left, Ti CPSL-27-7L (01)00841506143163  (01) 00841506143163	Clavicle Plate, Superior Lateral, 3.5mm, 8 Hole, Left, Ti CPSL-35-8L (01)00841506143187  (01) 00841506143187
Titanium Alloy Anterior Lateral Plates	
Clavicle Plate, Anterior Lateral, 2.7mm, 11 Hole, Ti CPAL-27-11 (01)00841506143194  (01) 00841506143194	Clavicle Plate, Anterior Lateral, 2.7/3.5mm, 12 Hole, Ti CPAL-2735-12 (01)00841506143200  (01) 00841506143200
Titanium Alloy Anterior Midshaft Plates	
Clavicle Plate, Anterior Midshaft, 2.7mm, 9 Hole, Ti CPAM-27-9 (01)00841506148946  (01) 00841506148946	Clavicle Plate, Anterior Midshaft, 3.5mm, 9 Hole, Ti CPAM-35-9 (01)00841506143217  (01) 00841506143217
Clavicle Plate Post	
Post, Clavicle Plate, Ti POST-CP (01)00841506149387  (01) 00841506149387	

Titanium Alloy PROTEAN Fragment Plates			
PROTEAN Fragment Plate, Anterior, 2.5mm, Ti PRT-CFS-ANT (01)00841506149226	 (01) 00841506149226	PROTEAN Fragment Plate, Curved Y, 2.5mm, Right, Ti PRT-CFS-CYR (01)00841506140568	 (01) 00841506140568
PROTEAN Fragment Plate, S, 2.5mm, Right, Ti PRT-CFS-SR (01)00841506140544	 (01) 00841506140544	PROTEAN Fragment Plate, Curved Y, 2.5mm, Left, Ti PRT-CFS-CYL (01)00841506140575	 (01) 00841506140575
PROTEAN Fragment Plate, S, 2.5mm, Left, Ti PRT-CFS-SL (01)00841506140551	 (01) 00841506140551		
2.7mm Titanium Alloy Fixed Angle Locking Screws			
Screw, Fixed Angle Locking, 2.7mm x 8mm, Ti FAL-27080-TS (01)00841506143347	 (01) 00841506143347	Screw, Fixed Angle Locking, 2.7mm x 17mm, Ti FAL-27170-TS (01)00841506144139	 (01) 00841506144139
Screw, Fixed Angle Locking, 2.7mm x 9mm, Ti FAL-27090-TS (01)00841506144092	 (01) 00841506144092	Screw, Fixed Angle Locking, 2.7mm x 18mm, Ti FAL-27180-TS (01)00841506143392	 (01) 00841506143392
Screw, Fixed Angle Locking, 2.7mm x 10mm, Ti FAL-27100-TS (01)00841506143354	 (01) 00841506143354	Screw, Fixed Angle Locking, 2.7mm x 19mm, Ti FAL-27190-TS (01)00841506144146	 (01) 00841506144146
Screw, Fixed Angle Locking, 2.7mm x 11mm, Ti FAL-27110-TS (01)00841506144108	 (01) 00841506144108	Screw, Fixed Angle Locking, 2.7mm x 20mm, Ti FAL-27200-TS (01)00841506143408	 (01) 00841506143408
Screw, Fixed Angle Locking, 2.7mm x 12mm, Ti FAL-27120-TS (01)00841506143361	 (01) 00841506143361	Screw, Fixed Angle Locking, 2.7mm x 22mm, Ti FAL-27220-TS (01)00841506143415	 (01) 00841506143415
Screw, Fixed Angle Locking, 2.7mm x 13mm, Ti FAL-27130-TS (01)00841506144115	 (01) 00841506144115	Screw, Fixed Angle Locking, 2.7mm x 24mm, Ti FAL-27240-TS (01)00841506143422	 (01) 00841506143422
Screw, Fixed Angle Locking, 2.7mm x 14mm, Ti FAL-27140-TS (01)00841506143378	 (01) 00841506143378	Screw, Fixed Angle Locking, 2.7mm x 26mm, Ti FAL-27260-TS (01)00841506143439	 (01) 00841506143439
Screw, Fixed Angle Locking, 2.7mm x 15mm, Ti FAL-27150-TS (01)00841506144122	 (01) 00841506144122	Screw, Fixed Angle Locking, 2.7mm x 28mm, Ti FAL-27280-TS (01)00841506143446	 (01) 00841506143446
Screw, Fixed Angle Locking, 2.7mm x 16mm, Ti FAL-27160-TS (01)00841506143385	 (01) 00841506143385	Screw, Fixed Angle Locking, 2.7mm x 30mm, Ti FAL-27300-TS (01)00841506143453	 (01) 00841506143453

Clavicle Fixation System Tray

Inventory Control Sheet

3.5mm Titanium Alloy Fixed Angle Locking Screws			
Screw, Fixed Angle Locking, 3.5mm x 8mm, Ti FAL-35080-TS (01)00841506145266	 (01) 00841506145266	Screw, Fixed Angle Locking, 3.5mm x 17mm, Ti FAL-35170-TS (01)00841506145358	 (01) 00841506145358
Screw, Fixed Angle Locking, 3.5mm x 9mm, Ti FAL-35090-TS (01)00841506145273	 (01) 00841506145273	Screw, Fixed Angle Locking, 3.5mm x 18mm, Ti FAL-35180-TS (01)00841506145365	 (01) 00841506145365
Screw, Fixed Angle Locking, 3.5mm x 10mm, Ti FAL-35100-TS (01)00841506145280	 (01) 00841506145280	Screw, Fixed Angle Locking, 3.5mm x 19mm, Ti FAL-35190-TS (01)00841506145372	 (01) 00841506145372
Screw, Fixed Angle Locking, 3.5mm x 11mm, Ti FAL-35110-TS (01)00841506145297	 (01) 00841506145297	Screw, Fixed Angle Locking, 3.5mm x 20mm, Ti FAL-35200-TS (01)00841506145389	 (01) 00841506145389
Screw, Fixed Angle Locking, 3.5mm x 12mm, Ti FAL-35120-TS (01)00841506145303	 (01) 00841506145303	Screw, Fixed Angle Locking, 3.5mm x 22mm, Ti FAL-35220-TS (01)00841506145402	 (01) 00841506145402
Screw, Fixed Angle Locking, 3.5mm x 13mm, Ti FAL-35130-TS (01)00841506145310	 (01) 00841506145310	Screw, Fixed Angle Locking, 3.5mm x 24mm, Ti FAL-35240-TS (01)00841506145426	 (01) 00841506145426
Screw, Fixed Angle Locking, 3.5mm x 14mm, Ti FAL-35140-TS (01)00841506145327	 (01) 00841506145327	Screw, Fixed Angle Locking, 3.5mm x 26mm, Ti FAL-35260-TS (01)00841506145440	 (01) 00841506145440
Screw, Fixed Angle Locking, 3.5mm x 15mm, Ti FAL-35150-TS (01)00841506145334	 (01) 00841506145334	Screw, Fixed Angle Locking, 3.5mm x 28mm, Ti FAL-35280-TS (01)00841506145464	 (01) 00841506145464
Screw, Fixed Angle Locking, 3.5mm x 16mm, Ti FAL-35160-TS (01)00841506145341	 (01) 00841506145341	Screw, Fixed Angle Locking, 3.5mm x 30mm, Ti FAL-35300-TS (01)00841506145488	 (01) 00841506145488
2.5mm Cobalt Chrome Alloy Variable Angle Locking Screws			
Screw, Variable Angle Locking, 2.5mm x 8mm, CoCr VAL-25080-CC (01)00841506136042	 (01) 00841506136042	Screw, Variable Angle Locking, 2.5mm x 17mm, CoCr VAL-25170-CC (01)00841506145600	 (01) 00841506145600
Screw, Variable Angle Locking, 2.5mm x 9mm, CoCr VAL-25090-CC (01)00841506136059	 (01) 00841506136059	Screw, Variable Angle Locking, 2.5mm x 18mm, CoCr VAL-25180-CC (01)00841506136134	 (01) 00841506136134
Screw, Variable Angle Locking, 2.5mm x 10mm, CoCr VAL-25100-CC (01)00841506136066	 (01) 00841506136066	Screw, Variable Angle Locking, 2.5mm x 19mm, CoCr VAL-25190-CC (01)00841506145617	 (01) 00841506145617
Screw, Variable Angle Locking, 2.5mm x 11mm, CoCr VAL-25110-CC (01)00841506136073	 (01) 00841506136073	Screw, Variable Angle Locking, 2.5mm x 20mm, CoCr VAL-25200-CC (01)00841506136141	 (01) 00841506136141

Clavicle Fixation System Tray

Inventory Control Sheet

Screw, Variable Angle Locking, 2.5mm x 12mm, CoCr VAL-25120-CC (01)00841506136080  (01) 00841506136080	Screw, Variable Angle Locking, 2.5mm x 22mm, CoCr VAL-25220-CC (01)00841506136158  (01) 00841506136158
Screw, Variable Angle Locking, 2.5mm x 13mm, CoCr VAL-25130-CC (01)00841506136097  (01) 00841506136097	Screw, Variable Angle Locking, 2.5mm x 24mm, CoCr VAL-25240-CC (01)00841506136165  (01) 00841506136165
Screw, Variable Angle Locking, 2.5mm x 14mm, CoCr VAL-25140-CC (01)00841506136103  (01) 00841506136103	Screw, Variable Angle Locking, 2.5mm x 26mm, CoCr VAL-25260-CC (01)00841506136172  (01) 00841506136172
Screw, Variable Angle Locking, 2.5mm x 15mm, CoCr VAL-25150-CC (01)00841506136110  (01) 00841506136110	Screw, Variable Angle Locking, 2.5mm x 28mm, CoCr VAL-25280-CC (01)00841506136189  (01) 00841506136189
Screw, Variable Angle Locking, 2.5mm x 16mm, CoCr VAL-25160-CC (01)00841506136127  (01) 00841506136127	Screw, Variable Angle Locking, 2.5mm x 30mm, CoCr VAL-25300-CC (01)00841506136196  (01) 00841506136196
2.7mm Cobalt Chrome Alloy Variable Angle Locking Screws	
Screw, Variable Angle Locking, 2.7mm x 8mm, CoCr VAL-27080-CC (01)00841506143460  (01) 00841506143460	Screw, Variable Angle Locking, 2.7mm x 17mm, CoCr VAL-27170-CC (01)00841506144245  (01) 00841506144245
Screw, Variable Angle Locking, 2.7mm x 9mm, CoCr VAL-27090-CC (01)00841506144207  (01) 00841506144207	Screw, Variable Angle Locking, 2.7mm x 18mm, CoCr VAL-27180-CC (01)00841506143514  (01) 00841506143514
Screw, Variable Angle Locking, 2.7mm x 10mm, CoCr VAL-27100-CC (01)00841506143477  (01) 00841506143477	Screw, Variable Angle Locking, 2.7mm x 19mm, CoCr VAL-27190-CC (01)00841506144252  (01) 00841506144252
Screw, Variable Angle Locking, 2.7mm x 11mm, CoCr VAL-27110-CC (01)00841506144214  (01) 00841506144214	Screw, Variable Angle Locking, 2.7mm x 20mm, CoCr VAL-27200-CC (01)00841506143521  (01) 00841506143521
Screw, Variable Angle Locking, 2.7mm x 12mm, CoCr VAL-27120-CC (01)00841506143484  (01) 00841506143484	Screw, Variable Angle Locking, 2.7mm x 22mm, CoCr VAL-27220-CC (01)00841506143538  (01) 00841506143538
Screw, Variable Angle Locking, 2.7mm x 13mm, CoCr VAL-27130-CC (01)00841506144221  (01) 00841506144221	Screw, Variable Angle Locking, 2.7mm x 24mm, CoCr VAL-27240-CC (01)00841506143545  (01) 00841506143545
Screw, Variable Angle Locking, 2.7mm x 14mm, CoCr VAL-27140-CC (01)00841506143491  (01) 00841506143491	Screw, Variable Angle Locking, 2.7mm x 26mm, CoCr VAL-27260-CC (01)00841506143552  (01) 00841506143552

Clavicle Fixation System Tray

Inventory Control Sheet

Screw, Variable Angle Locking, 2.7mm x 15mm, CoCr VAL-27150-CC (01)00841506144238  (01) 00841506144238	Screw, Variable Angle Locking, 2.7mm x 28mm, CoCr VAL-27280-CC (01)00841506143569  (01) 00841506143569
Screw, Variable Angle Locking, 2.7mm x 16mm, CoCr VAL-27160-CC (01)00841506143507  (01) 00841506143507	Screw, Variable Angle Locking, 2.7mm x 30mm, CoCr VAL-27300-CC (01)00841506143576  (01) 00841506143576
3.5mm Cobalt Chrome Alloy Variable Angle Locking-C Screws	
Screw, Variable Angle Locking-C, 3.5mm x 8mm, CoCr VALC-35080-CC (01)00841506143705  (01) 00841506143705	Screw, Variable Angle Locking-C, 3.5mm x 17mm, CoCr VALC-35170-CC (01)00841506144467  (01) 00841506144467
Screw, Variable Angle Locking-C, 3.5mm x 9mm, CoCr VALC-35090-CC (01)00841506144429  (01) 00841506144429	Screw, Variable Angle Locking-C, 3.5mm x 18mm, CoCr VALC-35180-CC (01)00841506143750  (01) 00841506143750
Screw, Variable Angle Locking-C, 3.5mm x 10mm, CoCr VALC-35100-CC (01)00841506143712  (01) 00841506143712	Screw, Variable Angle Locking-C, 3.5mm x 19mm, CoCr VALC-35190-CC (01)00841506144474  (01) 00841506144474
Screw, Variable Angle Locking-C, 3.5mm x 11mm, CoCr VALC-35110-CC (01)00841506144436  (01) 00841506144436	Screw, Variable Angle Locking-C, 3.5mm x 20mm, CoCr VALC-35200-CC (01)00841506145631  (01) 00841506145631
Screw, Variable Angle Locking-C, 3.5mm x 12mm, CoCr VALC-35120-CC (01)00841506143729  (01) 00841506143729	Screw, Variable Angle Locking-C, 3.5mm x 22mm, CoCr VALC-35220-CC (01)00841506145648  (01) 00841506145648
Screw, Variable Angle Locking-C, 3.5mm x 13mm, CoCr VALC-35130-CC (01)00841506144443  (01) 00841506144443	Screw, Variable Angle Locking-C, 3.5mm x 24mm, CoCr VALC-35240-CC (01)00841506145655  (01) 00841506145655
Screw, Variable Angle Locking-C, 3.5mm x 14mm, CoCr VALC-35140-CC (01)00841506143736  (01) 00841506143736	Screw, Variable Angle Locking-C, 3.5mm x 26mm, CoCr VALC-35260-CC (01)00841506145662  (01) 00841506145662
Screw, Variable Angle Locking-C, 3.5mm x 15mm, CoCr VALC-35150-CC (01)00841506144450  (01) 00841506144450	Screw, Variable Angle Locking-C, 3.5mm x 28mm, CoCr VALC-35280-CC (01)00841506145679  (01) 00841506145679
Screw, Variable Angle Locking-C, 3.5mm x 16mm, CoCr VALC-35160-CC (01)00841506143743  (01) 00841506143743	Screw, Variable Angle Locking-C, 3.5mm x 30mm, CoCr VALC-35300-CC (01)00841506145686  (01) 00841506145686
2.5mm Titanium Alloy Variable Angle Non-Locking Screws	
Screw, Non-Locking, 2.5mm x 8mm, Ti VAN-25080-TS (01)00841506137919  (01) 00841506137919	Screw, Non-Locking, 2.5mm x 17mm, Ti VAN-25170-TS (01)00841506145587  (01) 00841506145587
Screw, Non-Locking, 2.5mm x 9mm, Ti VAN-25090-TS (01)00841506137926  (01) 00841506137926	Screw, Non-Locking, 2.5mm x 18mm, Ti VAN-25180-TS (01)00841506138008  (01) 00841506138008

Clavicle Fixation System Tray

Inventory Control Sheet

Screw, Non-Locking, 2.5mm x 10mm, Ti VAN-25100-TS (01)00841506137933  (01) 00841506137933	Screw, Non-Locking, 2.5mm x 19mm, Ti VAN-25190-TS (01)00841506145594  (01) 00841506145594
Screw, Non-Locking, 2.5mm x 11mm, Ti VAN-25110-TS (01)00841506137940  (01) 00841506137940	Screw, Non-Locking, 2.5mm x 20mm, Ti VAN-25200-TS (01)00841506138015  (01) 00841506138015
Screw, Non-Locking, 2.5mm x 12mm, Ti VAN-25120-TS (01)00841506137957  (01) 00841506137957	Screw, Non-Locking, 2.5mm x 22mm, Ti VAN-25220-TS (01)00841506138022  (01) 00841506138022
Screw, Non-Locking, 2.5mm x 13mm, Ti VAN-25130-TS (01)00841506137964  (01) 00841506137964	Screw, Non-Locking, 2.5mm x 24mm, Ti VAN-25240-TS (01)00841506138039  (01) 00841506138039
Screw, Non-Locking, 2.5mm x 14mm, Ti VAN-25140-TS (01)00841506137971  (01) 00841506137971	Screw, Non-Locking, 2.5mm x 26mm, Ti VAN-25260-TS (01)00841506138046  (01) 00841506138046
Screw, Non-Locking, 2.5mm x 15mm, Ti VAN-25150-TS (01)00841506137988  (01) 00841506137988	Screw, Non-Locking, 2.5mm x 28mm, Ti VAN-25280-TS (01)00841506138053  (01) 00841506138053
Screw, Non-Locking, 2.5mm x 16mm, Ti VAN-25160-TS (01)00841506137995  (01) 00841506137995	Screw, Non-Locking, 2.5mm x 30mm, Ti VAN-25300-TS (01)00841506138060  (01) 00841506138060
2.7mm Titanium Alloy Variable Angle Non-Locking Screws	
Screw, Variable Angle Non-Locking, 2.7mm x 8mm, Ti VAN-27080-TS (01)00841506143224  (01) 00841506143224	Screw, Variable Angle Non-Locking, 2.7mm x 17mm, Ti VAN-27170-TS (01)00841506144023  (01) 00841506144023
Screw, Variable Angle Non-Locking, 2.7mm x 9mm, Ti VAN-27090-TS (01)00841506143989  (01) 00841506143989	Screw, Variable Angle Non-Locking, 2.7mm x 18mm, Ti VAN-27180-TS (01)00841506143279  (01) 00841506143279
Screw, Variable Angle Non-Locking, 2.7mm x 10mm, Ti VAN-27100-TS (01)00841506143231  (01) 00841506143231	Screw, Variable Angle Non-Locking, 2.7mm x 19mm, Ti VAN-27190-TS (01)00841506144030  (01) 00841506144030
Screw, Variable Angle Non-Locking, 2.7mm x 11mm, Ti VAN-27110-TS (01)00841506143996  (01) 00841506143996	Screw, Variable Angle Non-Locking, 2.7mm x 20mm, Ti VAN-27200-TS (01)00841506143286  (01) 00841506143286
Screw, Variable Angle Non-Locking, 2.7mm x 12mm, Ti VAN-27120-TS (01)00841506143248  (01) 00841506143248	Screw, Variable Angle Non-Locking, 2.7mm x 22mm, Ti VAN-27220-TS (01)00841506143293  (01) 00841506143293

Clavicle Fixation System Tray

Inventory Control Sheet

Screw, Variable Angle Non-Locking, 2.7mm x 13mm, Ti VAN-27130-TS (01)00841506144009  (01)00841506144009	Screw, Variable Angle Non-Locking, 2.7mm x 24mm, Ti VAN-27240-TS (01)00841506143309  (01)00841506143309
Screw, Variable Angle Non-Locking, 2.7mm x 14mm, Ti VAN-27140-TS (01)00841506143255  (01)00841506143255	Screw, Variable Angle Non-Locking, 2.7mm x 26mm, Ti VAN-27260-TS (01)00841506143316  (01)00841506143316
Screw, Variable Angle Non-Locking, 2.7mm x 15mm, Ti VAN-27150-TS (01)00841506144016  (01)00841506144016	Screw, Variable Angle Non-Locking, 2.7mm x 28mm, Ti VAN-27280-TS (01)00841506143323  (01)00841506143323
Screw, Variable Angle Non-Locking, 2.7mm x 16mm, Ti VAN-27160-TS (01)00841506143262  (01)00841506143262	Screw, Variable Angle Non-Locking, 2.7mm x 30mm, Ti VAN-27300-TS (01)00841506143330  (01)00841506143330
3.5mm Titanium Alloy Variable Angle Non-Locking Screws	
Screw, Variable Angle Non-Locking, 3.5mm x 8mm, Ti VAN-35080-TS (01)00841506143583  (01)00841506143583	Screw, Variable Angle Non-Locking, 3.5mm x 17mm, Ti VAN-35170-TS (01)00841506144351  (01)00841506144351
Screw, Variable Angle Non-Locking, 3.5mm x 9mm, Ti VAN-35090-TS (01)00841506144313  (01)00841506144313	Screw, Variable Angle Non-Locking, 3.5mm x 18mm, Ti VAN-35180-TS (01)00841506143637  (01)00841506143637
Screw, Variable Angle Non-Locking, 3.5mm x 10mm, Ti VAN-35100-TS (01)00841506143590  (01)00841506143590	Screw, Variable Angle Non-Locking, 3.5mm x 19mm, Ti VAN-35190-TS (01)00841506144368  (01)00841506144368
Screw, Variable Angle Non-Locking, 3.5mm x 11mm, Ti VAN-35110-TS (01)00841506144320  (01)00841506144320	Screw, Variable Angle Non-Locking, 3.5mm x 20mm, Ti VAN-35200-TS (01)00841506143644  (01)00841506143644
Screw, Variable Angle Non-Locking, 3.5mm x 12mm, Ti VAN-35120-TS (01)00841506143606  (01)00841506143606	Screw, Variable Angle Non-Locking, 3.5mm x 22mm, Ti VAN-35220-TS (01)00841506143651  (01)00841506143651
Screw, Variable Angle Non-Locking, 3.5mm x 13mm, Ti VAN-35130-TS (01)00841506144337  (01)00841506144337	Screw, Variable Angle Non-Locking, 3.5mm x 24mm, Ti VAN-35240-TS (01)00841506143668  (01)00841506143668
Screw, Variable Angle Non-Locking, 3.5mm x 14mm, Ti VAN-35140-TS (01)00841506143613  (01)00841506143613	Screw, Variable Angle Non-Locking, 3.5mm x 26mm, Ti VAN-35260-TS (01)00841506143675  (01)00841506143675
Screw, Variable Angle Non-Locking, 3.5mm x 15mm, Ti VAN-35150-TS (01)00841506144344  (01)00841506144344	Screw, Variable Angle Non-Locking, 3.5mm x 28mm, Ti VAN-35280-TS (01)00841506143682  (01)00841506143682
Screw, Variable Angle Non-Locking, 3.5mm x 16mm, Ti VAN-35160-TS (01)00841506143620  (01)00841506143620	Screw, Variable Angle Non-Locking, 3.5mm x 30mm, Ti VAN-35300-TS (01)00841506143699  (01)00841506143699

Clavicle Fixation System Tray

Inventory Control Sheet

2.5mm Titanium Alloy Big Head Variable Angle Non-Locking Screws			
Screw, Non-Locking, Big Head, 2.5mm x 8mm, Ti VANB-25080-TS (01)00841506143767	 (01)00841506143767	Screw, Non-Locking, Big Head, 2.5mm x 20mm, Ti VANB-25200-TS (01)00841506143828	 (01)00841506143828
Screw, Non-Locking, Big Head, 2.5mm x 10mm, Ti VANB-25100-TS (01)00841506143774	 (01)00841506143774	Screw, Non-Locking, Big Head, 2.5mm x 22mm, Ti VANB-25220-TS (01)00841506143835	 (01)00841506143835
Screw, Non-Locking, Big Head, 2.5mm x 12mm, Ti VANB-25120-TS (01)00841506143781	 (01)00841506143781	Screw, Non-Locking, Big Head, 2.5mm x 24mm, Ti VANB-25240-TS (01)00841506143842	 (01)00841506143842
Screw, Non-Locking, Big Head, 2.5mm x 14mm, Ti VANB-25140-TS (01)00841506143798	 (01)00841506143798	Screw, Non-Locking, Big Head, 2.5mm x 26mm, Ti VANB-25260-TS (01)00841506143859	 (01)00841506143859
Screw, Non-Locking, Big Head, 2.5mm x 16mm, Ti VANB-25160-TS (01)00841506143804	 (01)00841506143804	Screw, Non-Locking, Big Head, 2.5mm x 28mm, Ti VANB-25280-TS (01)00841506143866	 (01)00841506143866
Screw, Non-Locking, Big Head, 2.5mm x 18mm, Ti VANB-25180-TS (01)00841506143811	 (01)00841506143811	Screw, Non-Locking, Big Head, 2.5mm x 30mm, Ti VANB-25300-TS (01)00841506143873	 (01)00841506143873
Disposable Instruments			
K-Wire, 1.6 mm x 127mm KWIR-STD-15127 (01)00841506102504	 (01)00841506102504	K-Wire, 1.6 mm, Threaded, Olive KWIR-TOD-16127 (01)00841506141084	 (01)00841506141084
K-Wire, HCS, 1.6mm x 127mm, Double Trocar KWIR-DT-16127 (01)00841506109206	 (01)00841506109206		