
skeletal dynamics®

supinX™

DISTAL BICEPS REPAIR SYSTEM





Enabling anatomic reattachment of the distal biceps tendon


Anatomic by Design.
Purpose-Built for Repair.

The supinX™ Distal Biceps Repair System is designed to enable anatomic reattachment of the distal biceps tendon through a single-incision anterior approach using a minimally invasive, bone-preserving technique and proprietary automation technologies.

Anatomic

Facilitates biceps tendon reattachment to the native distal biceps footprint through either a single-head or double-head repair technique, supporting the distal biceps' natural biomechanics and function.







Anterior Single-Incision

Employs the standard and widely accepted single-incision anterior approach supporting a consistent and reproducible technique for accessing the radial tuberosity and reattachment of the biceps tendon.

Atraumatic

Utilizes a minimally invasive, 2mm bone-preserving tunnel to reattach the biceps tendon. The low-profile cortical button sits on the anterior aspect of the radial tuberosity, clear of neurovascular structures.





Automation

The patented Golden Retriever™ guide enables precise targeting and streamlined suture management. The proprietary Button Threader automatically loads the button onto the suture strands, simplifying tendon fixation.

Purpose-Built Instrumentation.

A optimized set of instruments designed to simplify the procedure and support a consistent, reproducible technique.



Golden Retriever™

Designed to facilitate passage and retrieval of the suture limbs through the 2mm bone tunnel, reducing soft tissue disruption and supporting a streamlined procedure.




Button Threader

Automates controlled delivery of the tendon suture limbs to the distal biceps button through the 2mm bone tunnel, helping to streamline button deployment and alleviate manual button to suture threading time.

Supports Natural Motion

By reattaching the biceps tendon to the native footprint, the supinX™ system helps recreate the distal biceps' natural biomechanics and function, including forearm supination strength and rotational range of motion.





INSTRUMENT SET INCLUDES	CATALOG #
Drill Guide	DBR-DG-20
Drill Guide Locking Rod, Distal Biceps Repair	DBR-DG-LR
Suture Passing Guide	DBR-SPG
STABLYX, Volar Capsule Elevator/Shoehorn	STX-VCE
Distal Biceps Drill, 2.0mm x 124mm	DRLL-20124
Instrument, Mini-Hohmann Retractor, Standard	INST-MHR-STD
Knot Pusher	KNT-PSH

STERILE IMPLANT KIT INCLUDES	CATALOG #
Distal Biceps Repair, Sterile Kit, Blue or **Black**	DBR-BLU/BLK-01
#2 Suture Loop, ST Needle	—
2-0 Suture Passing Needle Assembly	—
Button Threader Assembly	—
Circular Button, #2 Suture	—

Scan the QR Code
Learn more about our
supinX™ DISTAL BICEPS REPAIR SYSTEM



CE 2797
MKT-00510-00 V.2
October 2026

© 2026 Skeletal Dynamics
7300 North Kendall Drive, Suite 800, Miami, Florida 33156, 877-753-5396
Designed and Manufactured in the USA

EC REP
EMERGO EUROPE
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands